

Regular Cut Offs available for Drug Tests

Abbreviation	Description	Cut-off Level
AMP	Amphetamine	1000ng/mL
BAR	Barbiturates	300ng/mL
BUP	Buprenorphine	10ng/mL
BZO	Benzodiazepines	300ng/mL
COC	Cocaine	300ng/mL
mAMP	Methamphetamine	1000ng/mL
MDMA	MDMA Ecstasy	500ng/mL
MOP	Morphine	300ng/mL
MTD	Methadone	300ng/mL
OPI	Opiate	2000ng/mL
OXY	Oxycodone	100ng/mL
PCP	Phencyclidine	25ng/mL
PPX	Propoxyphene	300ng/mL
TCA	Tricyclic Antidepressants	1000ng/mL
THC	Marijuana/Cannabinoids	50ng/mL
TRA	Tramadol	200ng/mL
EDDP	2-ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine	300ng/mL
K2	JWH-018 / JWH-073	50ng/mL / 25ng/mL

Drug Test Available

CLIA WAIVED Test: AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OPI, OXY, PCP, TCA, THC

510K CLEARED And CE Marked Test: PPX

For Forensic Use only: TRA, EDDP, K2

Built in Adulteration parameters, offering visual check for urine integrity:

pH, Specific Gravity, Creatinine, Nitrite, Oxidant, Glutaraldehyde, Bleach, Pyridinium Chlorchromate

Wondfo USA Co., Ltd
545 Willowbrook Centre Parkway - Unit B
Willowbrook, IL 60527
Tel: 630-468-2199 Fax: 866-769-0168
Inquiry: wondfo@wondfousa.com
Website: www.wondfousa.com



**CLIA
WAIVED**