

Special Cut Offs available for Drug Tests

Abbreviation	Description	Cut-off Level
AMP	Amphetamine	500ng/mL
		300ng/mL
BZO	Benzodiazepines	200ng/mL
COC	Cocaine	150ng/mL
		100ng/mL
mAMP	Methamphetamine	500ng/mL
		300ng/mL
MOP	Morphine	100ng/mL
MTD	Methadone	200ng/mL
OXY	Oxycodone	300ng/mL
THC	Marijuana/Cannabinoids	20ng/mL
		40ng/mL
EDDP	2-ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine	100ng/mL
TRA	Tramadol	100ng/mL

Please inquire special order requirements for orders requesting Special Cut Offs

Drug Test Available

510K CLEARED And CE Marked Test:

AMP (300ng/mL), BZO(200ng/mL), COC(100ng/mL), mAMP(500ng/mL), mAMP(300ng/mL), MOP(100ng/mL), MTD(200ng/mL), THC(40ng/mL)

For Forensic Use only:

AMP (500ng/mL), COC(150ng/mL), OXY(300ng/mL), EDDP(100ng/mL), THC(20ng/mL), TRA(100ng/mL)

Built in Adulteration parameters, offering visual check for urine integrity:

pH, Specific Gravity, Creatinine, Nitrite, Oxidant, Glutaraldehyde, Bleach, Pyridinium Chlorchromate

