



# SAFE life™

## SAFElife T-Dip Multi-Drug Urine Test Panel

### Rapid Point-of-Care Toxicology Testing

#### Key Benefits:

- 2-16 configurations available
- Various formats and combinations
- Results in 5 minutes
- Customizable cut-off values (minimum order quantities required)
- Room temperature storage
- 24-month shelf life from date of manufacture
- Optional built-in adulteration for OEM products





















## CLIA-waived Panels

	Amphetamine (AMP)	1000, 500 ng/mL		Tricyclic Antidepressants (TCA)	1000 ng/mL
	Barbiturates (BAR)	300 ng/mL		Methamphetamine ( MET/mAMP)	1000, 500 ng/mL
	Benzodiazepines (BZO)	300 ng/mL		Methylenedioxymethamphetamine (MDMA)	500 ng/mL
	Buprenorphine (BUP)	10 ng/mL		Morphine (MOP/OPI300)	300 ng/mL
	Cocaine (COC)	300, 150 ng/mL		Opiate (OPI)	2000 ng/mL
	Cannabinoids (THC)	50 ng/mL		Oxycodone (OXY)	100 ng/mL
	2-ethylidene-1,5-dimethyl-3-3-diphenylpyrrolidine (EDDP)	300 ng/mL		Phencyclidine (PCP)	25 ng/mL
	Methadone (MTD)	300 ng/mL		Propoxyphene (PPX)	300 ng/mL



# SAFE life™

## Forensic Use Panels

	Monoacetylmorphine (6-MAM)	10 ng/mL		Synthetic Cannabisoids (K2)	50 ng/mL
	Amphetamine (AMP)	300, 500 ng/mL		Ketamine (KET)	300, 1000 ng/mL
	Benzodiazepines (BZO)	100, 200 ng/mL		Kratom (KRA)	300 ng/mL
	Cocaine (COC)	100, 150 ng/mL		Methamphetamine (MET/mAMP)	300, 500 ng/mL
	Cotinine (COT)	200 ng/mL		Morphine (MOP/OPI)	100 ng/mL
	EDDP	100 ng/mL		Methadone (MTD)	200 ng/mL
	Ethyl Glucuronide (EtG)	500 ng/mL		Methaqualone (MQL)	300 ng/mL
	Fentanyl (FTY)	20 ng/mL		Cannabinoids (THC)	25, 40 ng/mL
	Hydromorphone (HMO)	300 ng/mL		Tramadol (TRA)	100, 200, 1000 ng/mL

Generic stock adulteration testings include Creatinine (CR), Specific Gravity (S.G.) and pH. Other adulterates like Glutaraldehyde (GL), Oxidant (OX) & Nitrite (NI) are available through special order.



# SAFE life™

## T-Dip Single Drug Urine Test Card

⌘ Kit Size: 25 / Case Size: 20

Item Code	Description	Complexity
WDBU-114	1-Panel Card; BUP Buprenorphine	CLIA WAIVED
WDBZ-114	1-Panel Card; BZO Benzodiazepines	CLIA WAIVED
WDCC-114	1-Panel Card; COC Cocaine	CLIA WAIVED
WDMA-114	1-Panel Card; mAMP/MET Methamphetamine	CLIA WAIVED
WDMT-114	1-Panel Card; MTD Methadone	CLIA WAIVED
WDOP-114	1-Panel Card; OPI Opiate 2000	CLIA WAIVED
WDOX-114	1-Panel Card; OXY Oxycodone	CLIA WAIVED
WDPC-114	1-Panel Card; PCP Phencyclidine	CLIA WAIVED
WDTH-114	1-Panel Card; THC Marijuana/Cannabinoids	CLIA WAIVED
WCOT-114	1-Panel Card; COT Cotinine	Forensic Use Only
WDOA-K2	1-Panel Card; K2 Synthetic Cannabinoid	Forensic Use Only
WDTR-114	1-Panel Card; TRA Tramadol	Forensic Use Only
WETG-114	1-Panel Card; EtG Ethyl Glucuronide	Forensic Use Only
WETG-300	1-Panel Card; EtG Ethyl Glucuronide	Forensic Use Only
WFTY-114	1-Panel Card; FTY Fentanyl	Forensic Use Only
WKRA-114	1-Panel Card; KRA Kratom	Forensic Use Only
W6AM-114	1-Panel Card; 6-MAM Monoacetylmorphine	Forensic Use Only
WMQL-114	1-Panel Card; MQL Methaqualone	Forensic Use Only



# SAFE life™

## T-Dip Multi-Drug Urine Test Panel

⌘ Kit Size: 25 / Case Size: 20

Item Code	Description	Complexity
WDOA-234	3-Panel Card; COC, mAMP, THC	CLIA WAIVED
WDOA-154	5-Panel Card; AMP, COC, OPI, PCP, THC	CLIA WAIVED
WDOA-254	5-Panel Card; AMP, COC, mAMP, OPI, THC	CLIA WAIVED
WDOA-354	5-Panel Card; COC, mAMP, OPI, OXY, THC	CLIA WAIVED
WDOA-554	5-Panel Card; COC, mAMP, MDMA, OPI, THC	CLIA WAIVED
WDOA-654	5-Panel Card; BZO, COC, mAMP, OPI, THC	CLIA WAIVED
WDOA-754	5-Panel Card; AMP, BZO, COC, OPI, THC	CLIA WAIVED
WDOA-264	6-Panel Card; AMP, BZO, COC, mAMP, OPI, THC	CLIA WAIVED
WDOA-564	6-Panel Card; BZO, COC, mAMP, OPI, OXY, THC	CLIA WAIVED
WDOA-865	6-Panel Card; BUP, BZO, COC, MOP, MTD, OXY	CLIA WAIVED
WDOA-274	7-Panel Card; AMP, BZO, COC, mAMP, OPI, OXY, THC	CLIA WAIVED
WDOA-295	9-Panel Card; AMP, BUP, BZO, COC, MOP, MTD, OXY, PPX, THC	CLIA WAIVED
WDOA-994NT	9-Panel Card without THC; AMP, BAR, BUP, BZO, COC, MET, MTD, OPI, OXY	CLIA WAIVED
WDOA-1104	10-Panel Card; AMP, BAR, BZO, COC, mAMP, MTD, OPI, PCP, TCA, THC	CLIA WAIVED
WDOA-3104	10-Panel Card; AMP, BAR, BZO, COC, mAMP, MDMA, MTD, OPI, PCP, THC	CLIA WAIVED
WDOA-4104	10-Panel Card; AMP, BAR, BZO, COC, mAMP, MTD, OPI, OXY, PCP, THC	CLIA WAIVED
WDOA-7104	10-Panel Card; BAR, BUP, BZO, COC, mAMP, MTD, OPI, OXY, TCA, THC	CLIA WAIVED
WDOA-9104	10-Panel Card; AMP, BAR, BUP, BZO, COC, mAMP, MTD, OPI, OXY, THC	CLIA WAIVED
WDOA-1124	12-Panel Card; AMP, BAR, BZO, COC, mAMP, MDMA, MTD, OPI, OXY, PCP, PPX, THC	CLIA WAIVED
WDOA-3124	12-Panel Card; AMP, BAR, BUP, BZO, COC, MDMA, MTD, TCA, OPI, OXY, PCP, THC	CLIA WAIVED
WDOA-6124	12-Panel Card; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MTD, OPI, OXY, PCP, THC	CLIA WAIVED
WDOA-6125	12-Panel Card; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, THC	CLIA WAIVED
WDOA-7125	12-Panel Card; AMP, BAR, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, TCA, THC	CLIA WAIVED
WDOA-8145EF	14-Panel Card; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, THC, ETG, FTY	Forensic Use Only