



SAFE life™

C-Cup Multi-Drug Urine Test Cup

Rapid Point-of-Care Toxicology Testing

- Cost Effective
- Reduce storage space
- Lower shipping expense
- CLIA Waived, 510(k) cleared & Forensic Use Only
- Same technical features as the T-Cup
- Test up to 20 drugs per cup



CLIA-waived Drugs & Cutoff Levels

AMP AMP AMP	Amphetamine (AMP)	1000, 500 ng/mL	TCA TCA TCA	Tricyclic Antidepressants (TCA)	1000 ng/mL
BAR BAR BAR	Barbiturates (BAR)	300 ng/mL	MET MET MET	Methamphetamine (MET/mAMP)	1000, 500 ng/mL
BZO BZO BZO	Benzodiazepines (BZO)	300 ng/mL	MDMA MDMA MDMA	Methylenedioxymethamphetamine (MDMA)	500 ng/mL
BUP BUP BUP	Buprenorphine (BUP)	10 ng/mL	MOP MOP MOP	Morphine (MOP/OPI300)	300 ng/mL
COC COC COC	Cocaine (COC)	300, 150 ng/mL	OPI OPI OPI	Opiate (OPI)	2000 ng/mL
THC THC THC	Cannabinoids (THC)	50 ng/mL	OXY OXY OXY	Oxycodone (OXY)	100 ng/mL
EDDP EDDP EDDP	2-ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine (EDDP)	300 ng/mL	PCP PCP PCP	Phencyclidine (PCP)	25 ng/mL
MTD MTD MTD	Methadone (MTD)	300 ng/mL	PPX PPX PPX	Propoxyphene (PPX)	300 ng/mL

Forensic Use Panels

6-MAM 6-MAM 6-MAM	Monoacetylmorphine (6-MAM)	10 ng/mL	K2 K2 K2	Synthetic Cannabinoid (K2)	50 ng/mL
AMP AMP AMP	Amphetamine (AMP)	300, 500 ng/mL	KET KET KET	Ketamine (KET)	300, 1000 ng/mL
BUP BUP BUP	Buprenorphine (BUP)	5 ng/mL	KRA KRA KRA	Kratom (KRA)	100, 300 ng/mL
BZO BZO BZO	Benzodiazepines (BZO)	100, 200 ng/mL	LSD LSD LSD	Lysergic Acid Diethylamide (LSD)	20 ng/mL
COC COC COC	Cocaine (COC)	100, 150 ng/mL	MET MET MET	Methamphetamine (MET/mAMP)	300, 500 ng/mL
COT COT COT	Cotinine (COT)	200 ng/mL	MOP MOP MOP	Morphine (MOP or OPI)	100 ng/mL
EDDP EDDP EDDP	EDDP	100 ng/mL	MTD MTD MTD	Methadone (MTD)	200 ng/mL
EtG EtG EtG	Ethyl Glucuronide (EtG)	500 ng/mL	MQL MQL MQL	Methaqualone (MQL)	300 ng/mL
FTY FTY FTY	Fentanyl (FTY)	20 ng/mL	PGB PGB PGB	Pregabalin (PGB)	500 ng/mL
GAB GAB GAB	Gabapentin (GAB)	2000 ng/mL	THC THC THC	Cannabinoids (THC)	15, 25, 40 ng/mL
HMO HMO HMO	Hydromorphone (HMO)	300 ng/mL	TRA TRA TRA	Tramadol (TRA)	100, 200, 1000 ng/mL

Generic stock adulteration testings include Creatinine (CR), Specific Gravity (S.G.) and pH. Other adulterates like Glutaraldehyde (GL), Oxidant (OX) & Nitrite (NI) are available through special order.



SAFE life™

Product Information

※ Kit Size: 25 / Case Size: 4

Item Code	Description	Complexity
CDOA-254	5-Panel C-Cup; AMP, COC, mAMP, OPI, THC	CLIA WAIVED
CDOA-564	6-Panel C-Cup; BZO, COC, mAMP, OPI, OXY, THC	CLIA WAIVED
CDOA-8105	10-Panel C-Cup; AMP, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, THC	CLIA WAIVED
CDOA-3124	12-Panel C-Cup; AMP, BAR, BUP, BZO, COC, MDMA, MTD, OPI, OXY, PCP, TCA, THC	CLIA WAIVED
CDOA-6125	12-Panel C-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, THC	CLIA WAIVED
CDOA-6125A3	12-Panel C-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, THC with 3AD	CLIA WAIVED
CDOA-L6125	12-Panel C-Cup, low cutoff; AMP500, BAR, BUP, BZO, COC150, mAMP500, MDMA, MOP, MTD, OXY, PCP, THC	CLIA WAIVED
CDOA-9125A3E3F	12-Panel C-Cup; AMP, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, THC, ETG300, FTY with 3AD	Forensic Use Only
CDOA-L9135NTA3	13-Panel C-Cup without THC, low cutoff; AMP500, BAR, BUP, BZO, COC150, MDMA, MET500, MOP, MTD, OXY, PCP, PPX, TCA with 3AD	CLIA WAIVED
CDOA-9165EFTK	16-Panel C-Cup; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, THC, ETG, FTY, TRA, K2	Forensic Use Only
CDOA-9165EFTKA3	16-Panel C-Cup; AMP, BAR, BUP, BZO, COC, MDMA, MET, MOP, MTD, OXY, PCP, THC, ETG, FTY, TRA, K2 with 3AD	Forensic Use Only
CDOA-S1185EFTKK	18-Panel C-Cup, Special TRA; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, TCA, THC, ETG, FTY, TRA100, K2, KRA	Forensic Use Only

Rel: 2024/05/17